

# Chemical and Puncture Resistant gloves

1. Remove exam gloves, perform hand-hygiene and dismiss Patient
2. Go back to operatory and put on utility gloves – hanging in operatory
3. Dispose of all trash and sharps and get instruments ready to transport
4. Transport instruments to sterilization area in closed, locked container and deposit them in ultrasonic, instrument washer or holding area (spray with enzymatic gel)
5. Perform any tasks in sterilization area with Utility Gloves on
6. Wash Utility Glove with soap and water and dry with paper towel or wipe with surface disinfectant
7. Return to operatory and properly clean surfaces
8. Use TB Killing disinfectant to wipe down utility gloves before they are removed.
9. Hang gloves to dry

